

# HORIZON

## HORIZON 6, 12, & 24 FLEX

### **PROGRAMMABLE ROUTINE CENTRIFUGES**

Satisfy the most diverse processing requirements with HORIZON Flex routine centrifuges. Customize settings and check every detail on the digital display, or pre-program up to 10 cycles.

#### **Easy to Use**

The fully digital display makes it easy to customize cycle parameters on the fly, or set and select one of the 10 programmable cycles. Spin 75-100 mm (3-10 mL) tubes without switching tube holders or adding cushions.

#### **Easy to Monitor**

Process samples more efficiently with HORIZON's LED lid lighting: off when the centrifuge is ready to load, on when running, and flashing at the end of cycle. Once the unit is started, the display shows remaining spin time.

#### **Built to Last**

No routine maintenance. Brushless motor, engineered composite components.

#### **Engineered for Safety**

Protect your lab with imbalance detection and a shatterproof lid which can only be opened when the rotor is stopped, even if the unit loses power.

#### **A Great Value**

Unbeatable quality for any budget. Each HORIZON centrifuge is ready to run with horizontal rotor, universal tube holders, and 2-year warranty included in the affordable price.



Tube Capacity	6, 12, or 24 (75-100 mm/3-10 mL)
Applications	Chemistry, Coag, PPP, PRP, Urine
Centrifugation	Horizontal (max 2,000 xg/3,800 RPM)
Included	Rotor, buckets and tube holders
Warranty *	2 years (rotor: lifetime warranty)

\* *Extended warranty available*



**HORIZON 6 Flex**



**HORIZON 12 Flex**



**HORIZON 24 Flex**

**Designed, built, and supported in the USA**

# HORIZON 6, 12, & 24 FLEX



## Free 30-Day Evaluation (U.S. Only)

Part Numbers	HORIZON 6 Flex	HORIZON 12 Flex	HORIZON 24 Flex
Centrifuge	00-376-009-000	00-383-009-000	00-384-009-000
Replacement 75-100 mm tube holder/Bucket	7713079 (1 qty)	7713079 (1 qty)	7713023 (1 qty)
50 mL tube holder			7713037
Cap			7713035
SmartView Platform	00-079-009-001		
1 yr Warranty Extension	03-0-0007-0001	03-0-0007-0001	03-0-0007-0001
2 yr Warranty Extension	03-0-0007-0002	03-0-0007-0002	03-0-0007-0002
3 yr Warranty Extension	03-0-0007-0006	03-0-0007-0006	03-0-0007-0006
Custom Settings *	00-100-100-101	00-100-100-101	00-100-100-101

\* Add this part number to your centrifuge order to have your custom validated cycles pre-programmed at the factory

Specifications	HORIZON 6 Flex	HORIZON 12 Flex	HORIZON 24 Flex
Sound Level	64 dBA	56 dBA	59 dBA
Max G-Force/Speed	2,000 xg @ 3,800 RPM	2,000 xg @ 3,700 RPM	2,000 xg @ 3,400 RPM
Dimensions (W x D x H)	12 x 14 x 9 in 30 x 36 x 23 cm	13 x 15 x 9 in 33 x 38 x 23 cm	15 x 17 x 9 in 38 x 43 x 23 cm
Weight	12 lbs (5.4 kg)	30 lbs (13.7 kg)	39 lbs (17.7 kg)
Power	95 to 253 V – 50/60 Hz		

FDA LISTED



**RoHS**  
Compliant



**ISO13485**  
certified



druckerdiagnostics.com  
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)  
+1-814-692-7661

